

Ocutears® Alo+

Instructions for use
Preservative free ophthalmic solution

What is Ocutears® Alo+ and what it is used for?

Ocutears® Alo+ is based on cross-linked sodium hyaluronate (0.2%) and aloe vera (0.1%), formulated for the treatment of dry eye symptoms, ocular discomfort and corneal protection.

Cross-linked sodium hyaluronate gives the product a high viscoelasticity ensuring a high degree of hydration and permanence on the ocular surface, providing long lasting relief from dry eye symptoms.

Aloe Vera is known for moisturising, decongestant and lubricating properties, giving the solution a higher moisturising capacity, therefore improving visual comfort and prolonging the feeling of freshness. Aloe Vera has also proven to be useful in case of redness or particularly sensitive eyes, with highly moisturising and long-lasting protective action.

Ocutears® Alo+ is indicated

- for the relief of dry eye symptoms, irritation and redness caused by environmental conditions or medications.
- can also be administered for ocular lubrication and comfort when using contact lenses

Ocutears® Alo+ is useful during and after:

- prolonged exposure to sunlight;
- Exposure to air-conditioning;
- prolonged use of digital screens;


Ocutears® Alo+ could also relieve the post surgery discomfort of the eye

How to use Ocutears® Alo+

Administer 1 or 2 drops in each eye one or more times a day or as advised by your doctor or pharmacist. Wash your hands before opening the container and pay attention to hygiene during handling.

Bottle:


- 1) Before first opening pull and remove the tamper evident band completely.
- 2) Gently squeeze 1 drop into the eye.
- 3) Close firmly the bottle after use.



1




2



3

UD:

- 1) Take the strip of mini containers and remove one.
- 2) Open the container by twisting the cap.
- 3) Gently squeeze 1 drop into the eye.
- 4) Discard the remaining product in the container after use.



Read the whole leaflet carefully before you start using this product.
 Keep this leaflet; you may need to read it again.
 If you have any further questions, ask your doctor or pharmacist.
 If symptoms persist or you have any side effects please consult your doctor.

Warnings:

- Do not use the product if allergic or hypersensitive to any of the components;
- **In the case of serious corneal pathologies, consult a doctor before use;**
- Use the bottle within 90 days from opening; the single dose container should be used immediately after opening and discarded after use;
- Do not use the product after the expiry date indicated on the pack;
- Do not use the product if the packaging is damaged;
- Keep away from heat, sunlight, moisture and the reach of children;
- Do not ingest;
- Take care that the tip of the container does not come into contact with the eye;
- Store the product below 25°C. Do not freeze;


Ocutears® Alo+ contains:

Cross-linked sodium hyaluronate 0.2%, Aloe Vera gel 0.1%, Disodium tetraborate decahydrate, Boric acid, Sodium chloride, Water for injections (WFI).

Presentations: 10 ml bottle (PFMD), 0.35 ml x 15 UD, 0.35 ml x 5 UD (free sample)

For further information please contact your doctor or your local Santen Office at www.santen.eu

 OFFHEALTH S.p.A., Via Giovanni Paisiello 10, 50144 Florence, Italy



This leaflet was revised in **March 2021**



										
Manufacturing date	Manufacturer	Store below 25°C	Keep away from sunlight	Sterilized using aseptic processing techniques	Batch code	Use by date	Consult instructions for use	Do not re-use	Do not use if package is damaged	Keep dry

carton OSD

Ocutears®Alo+

for irritated Dry Eye
Preservative free artificial tear
HA 0.2% & Aloe Vera Highly moisturising
10 ml

Lubricating ophthalmic solution indicated for the relief of dry eye symptoms, irritation and redness caused by environmental conditions or medications.

Composition: Cross-linked sodium hyaluronate 0.2%, Aloe Vera gel 0.1%, Disodium tetraborate decahydrate, Boric acid, Sodium chloride, WFI.

 OFFHEALTH S.p.A., Via Giovanni Paisiello 10, 50144 Florence, Italy



Santen    + EAN13 **REF** 0560H

EYE HEALTH SPECIALISTS

Rev. 07 of 12/2020

label OSD

Ocutears®Alo+

Preservative free artificial tear
HA 0.2%& Aloe Vera, Highly moisturising
10 ml

 OFFHEALTH S.p.A., Florence, Italy



Santen   

carton UD SALE

Ocutears®Alo+

for irritated Dry Eye
HA 0.2% & Aloe Vera

Highly moisturising preservative free artificial tear

Lubricating ophthalmic solution indicated for the relief of dry eye symptoms, irritation and redness caused by environmental conditions or medications.

Composition: Cross-linked sodium hyaluronate 0.2%, Aloe Vera gel 0.1%, Disodium tetraborate decahydrate, Boric acid, Sodium chloride, WFI.

 OFFHEALTH S.p.A., Via Giovanni Paisiello 10, 50144 Florence, Italy


1370

STERILE A



Santen

LOT



+ EAN13

REF 0500H

15
Single-dose
Containers
of 0.35 ml
Preservative free

EYE HEALTH SPECIALISTS

Rev. 07 of 12/2020

carton UD SAMPLE

Ocutears®Alo+

for irritated Dry Eye
HA 0.2% & Aloe Vera

Highly moisturising preservative free artificial tear

Lubricating ophthalmic solution indicated for the relief of dry eye symptoms, irritation and redness caused by environmental conditions or medications.

Composition: Cross-linked sodium hyaluronate 0.2%, Aloe Vera gel 0.1%, Disodium tetraborate decahydrate, Boric acid, Sodium chloride, WFI.

 OFFHEALTH S.p.A., Via Giovanni Paisiello 10, 50144 Florence, Italy



STERILE A



Santen

LOT



+ EAN13

REF

051OH

5
Single-dose
Containers
of 0.35 ml
Preservative free

EYE HEALTH SPECIALISTS
FREE SAMPLE NOT FOR SALE

Rev. 07 of 12/2020